


Bacti-cinerator Sterilizer KT-IS800A



Description:

KT-IS800A Bacti-cinerator Sterilizer is used in degassing for inoculating loops and needles, instead of alcohol lamps. Due to its convenience, sterilizer is widely used in the environment of biosafety cabinets, clean benches. They are also used besides exhaust fans or on mobile vans.

Features:

- ◆ Safety, convenient. Totally take the place of alcohol lamps.
- ◆ Only 5 seconds to 7 seconds are needed, because of its highest temperature of 825°C.
- ◆ Smart, light, easy to handle.
- ◆ Used in anaerobic chamber is allowed.
- ◆ Podzolic organics prevent pollutions effectively.
- ◆ Smart technique makes it possible to use for a long time and much safer.

Application:

Widely used in biosafety cabinets, clean bench, near exhaust fan, mobile vehicles, even outdoor and other severe environment can be easily used.

Parameters:

Parameter	Type	KT-IS800A
Max. temp. in central zone		825°C±25°C
Max. outer diameter of sterilization goods		Φ15mm
Length of heating zone		140mm
Heating time		15min
Ambient temperature		5°C~40°C
Fuse		250V 3A Φ5×20
Voltage		AC220V or AC110V, 50-60Hz
Power		150W
Dimension(mm)		120×155×180
Weight(kg)		1.0